## CS 530

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### Introduction

Dear Client,

Thank you for purchasing the Pic CS 530 automatic wrist blood pressure monitor.

This fully automatic monitor is easy to use and ideal for daily measurements. The large display clearly shows systolic and diastolic pressure and pulse rate readings, as well as other useful information.

With Pic CS 530 you can store up to 60 measured values for two different users, making it ideal for people who want to regularly monitor their blood pressure readings.

The CS 530 is compact and portable, making it perfect for home use and travel.

Blood pressure measurements taken with this monitor are equivalent to those obtained by a trained person using the cuff/stethoscope auscultation method, within the limits prescribed by the American National Standard for Electronic or Automated Sphygmomanometers.

The monitor is designed for home use by adults.

 Note: Please read this manual carefully before use and keep for future reference. For more information about your own blood pressure, please consult your doctor.

### Features.:

- Smart Sense technology
- One Touch easy measurement
- pressure level indicator
- irregular heart beat indicator
- customisable dual memory (60 + 60 readings can be stored)
- average of last three readings
- date and time
- large, easy to read display


### Quick guide

A. insert the batteries  
Type of batteries: 2 alkaline LR03  
AAA batteries

B. applying the cuff  
Wrap the cuff around your left  
wrist, with the palm of your hand  
facing upwards

C. correct position  
Rest your wrist on a surface, level  
with your heart

**D. taking a reading**  
Press the ON/OFF button and keep still during the measurement: the slightest movement could result in an inaccurate measurement

**E. automatic memory**  
The symbol "mem." will flash after the reading has been taken and stored

**F. displaying stored data**  
Press the Select User button, then the M button to retrieve previously stored readings

For more details on each operation, please read the relative sections in this manual.

### **Preliminary remarks**

This blood pressure monitor conforms to European standards and bears the CE mark "CE 0068".

The quality of the device has been verified and conforms to the provisions of the European Council directive 93/42/EEC of 14 June 1993 on medical devices, as well as the directive 89/336/EEC on electromagnetic compatibility:

**EN 1060-1: 1995/A1: 2002** Non-invasive blood pressure measuring equipment. General requirements.

**EN 1060-3 (Feb. 1997)** Non-invasive blood pressure measuring equipment. Supplementary requirements for electromechanical blood pressure measuring systems.

To ensure accurate measurements, the calibration of the unit should be checked every two years.

### **How the Pic CS 530 works**

The Pic CS 530 blood pressure monitor uses the oscillometric method to measure arterial blood pressure. Before the cuff is inflated, the device will establish a baseline cuff pressure equivalent to the air pressure.

At the press of the button, the cuff is automatically inflated.

As it inflates, the monitor will take a first arterial pressure reading, analysing the oscillations.

Once a pressure of 180 mmHg has been reached, the cuff will start to slowly deflate, accurately measuring the amplitude and gradient of the pressure oscillations; once the cuff has entirely deflated, the systolic and diastolic pressure and pulse rate are displayed.

### **Arterial blood pressure**

*What is arterial blood pressure?*

Arterial blood pressure is the pressure exerted on the artery as blood flows through it.

The pressure measured when the heart contracts and pumps out blood is called systolic (or highest) blood pressure, while the pressure measured when the heart dilates and blood flows back is called diastolic (lowest) blood pressure.

### Why measure your blood pressure?


High blood pressure is by far the most common health problem afflicting modern man and can lead to a number of complications. Taking your own blood pressure is therefore a useful means of prevention and of controlling your health. You should however consult your doctor, who will be able to tell you what your ideal blood pressure should be.

### Blood pressure fluctuation

Blood pressure fluctuates all the time!

Blood pressure varies during the month and even throughout the day. It is also influenced by the season and temperature.

The graph below shows, by way of example, normal fluctuations that may occur during the day.

### Names/function of the various parts:

data transfer port reserved for Technical Service Centre. Warning: do not insert any object, doing so could interfere with the accuracy of measurements.

Date/time adjustment buttons

LCD display

Battery cover

Memory button


2 AAA 1.5V batteries

Cuff

User selection button

ON/OFF button

## Information shown on the display

## Meaning of the symbols

	User 1	appears when the device is used by User 1
	User 2	appears when the device is used by User 2
	Flat battery	appears when the battery needs to be replaced
	"Memory"	indicates the consecutive reading stored in the monitor
	"Average"	indicates the average value of the last few readings (maximum 3 readings)
	Pressure level indicator	compares the pressure value measured with values indicated by the Joint National Committee on Prevention, Detection, Evaluation and Treatment of High Blood Pressure (see below)
	Irregular heart beat indicator	informs the user that the heart beat is irregular during the measurement

### Pressure level indicator based on the classification of the Joint National Committee (USA)

The 7th report of the Joint National Committee on Prevention, Detection, Evaluation and Treatment of High Blood Pressure (J.N.C.) has identified 4 arterial blood pressure groups which indicate high blood pressure risk values.

The Pic CS 530 monitor has an innovative indicator which classifies, and displays, blood pressure according to the 4 groups outlined by the J.N.C. and gives an initial indication of the presumed risk (pre-hypertension, level 1 hypertension, level 2 hypertension)

Systolic pressure	<120 and Diastolic pressure <80	120-139 or 80-89	140-159 or 90-99	>160 or >100
J.N.C. classification	Normal	Pre-hypertension	Level 1 hypertension	Level 2 hypertension
Symbol on the display	no symbol			

You should however consult your doctor, he or she will be able to tell you your normal blood pressure range as well as the point at which you should consider yourself at risk.

The J.N.C. classification is only a general guideline, as individual blood pressure measurements vary according to population, gender, age group, etc.

### Irregular heart beat indicator

The monitor has an irregular heart beat indicator.

When taking readings, the monitor calculates the heart beat frequency and if it detects a beat or beats that are far higher or lower than the mean frequency, a symbol appears on the display. 

If this symbol appears regularly, we suggest you to contact your doctor. An irregular heart beat may be due to numerous factors and only your doctor can assess the situation properly.

The indicators "Pressure level" and "Irregular heart beat" have no diagnostic value but simply provide general information.

### Power supply

Inserting the batteries

1. push and slide off the battery cover in the direction of the arrow.  
Batteries: 2 alkaline LRO3 AAA batteries

2. insert the first battery in the lower position on top of the ribbon.

3. insert the other battery and fold the ribbon over it.

4. replace the cover; make sure it closes with a click.

Batteries should be replaced when:


- the flat battery symbol appears on the display
- nothing appears on the display when you press the ON/OFF button

Always replace both batteries at the same time.


Remove the batteries from the monitor if it is not used for long periods.

### Applying the monitor


Remove all watches, jewellery, etc. before applying the cuff to bare skin for a correct measurement.


**1.** wrap the cuff around your left wrist, with the palm facing upwards as shown in figure A.

**2.** make sure the edge of the cuff is about 1 cm from your palm as shown in figure B.


**3.** in order to ensure accurate measurements, fasten the Velcro strap securely around the wrist, as shown in figure C.

**4.** If a doctor has diagnosed poor circulation in your left arm, place the cuff around your right wrist, as shown in figure D.

### Correct position for measurements

**1.** rest you elbow on a table so the cuff is level with your heart, as shown in figure A.

Note: Your heart is located slightly below your armpit. Relax your entire body, especially the area between your elbow and fingers.

Fig. B

2. if the cuff is not at the same level as your heart or if you are unable to keep your arm completely still during the measurement, use a soft object such as a folded towel, to support you arm, as shown in figure B. Do not allow hard objects to come into contact with the cuff.

3. turn your palm upwards.

4. sit upright and take 5 to 6 deep breaths before taking the reading.

Fig. C

Do not lean back while the measurement is being taken, as shown in figure C.

### Tips for easier readings


- do not measure your blood pressure immediately after eating a large meal. To obtain more accurate readings, wait at least one hour before measuring your blood pressure.
- do not smoke or consume alcohol before measuring your blood pressure.
- do not measure your blood pressure if you are tired.
- it is important to relax during the measurement. Try to rest for 15 minutes before measuring your blood pressure.
- do not talk or move your arm or hand muscles when measuring your blood pressure.
- measure your blood pressure at normal body temperature; if you feel hot or cold, wait a while before taking a measurement.
- if the monitor is stored at a very low temperature (near freezing), put it in a warm place for at least an hour before use
- leave at least 5 minutes between each measurement.

### Taking a measurement

1. press the Select User button to choose user **1** or user **2**

2. wrap the cuff around your wrist and press the **ON/OFF** button: all the digits will light up to check the display functions. This will take about 2 seconds.


After all symbols have appeared on the display, a flashing 0 will appear: the monitor is now ready to take a measurement.

3. once the monitor has inflated the cuff to a pressure of approximately 180 mmHg, the cuff will slowly deflate and the measurement will be taken. If the cuff reaches a pressure of 300 mmHg when inflating, the monitor will immediately start to deflate it for safety reasons.

4. after the measurement has been taken, the systolic and diastolic pressure and pulse rate values will be shown at the same time on the display and automatically saved in the previously selected user memory.

5. press the ON/OFF button to turn off the monitor. The monitor will automatically switch off after 60 seconds, if no buttons are pressed.

#### Notes

a) the monitor will automatically reinflate the cuff to a pressure of approximately 200 mmHg and take another reading, if it detects the user needs greater pressure to measure blood pressure.

b) press the ON/OFF button or the M button to stop a reading: the cuff will immediately deflate

#### Retrieving stored values

1. the monitor has two User Memories. Each one can store up to 60 readings.


2. to see stored values, press the Select User Memory button to choose the User Memory (1 or 2) press the M button to display stored values

3. the first digit to appear is the average value of the last few readings (a maximum of three); press the M button again to display previously stored readings, starting with the most recent.


## Deleting stored values

1. to delete stored values, press the Select User Memory button to choose the User Memory (1 or 2)
2. press down the M button for 5 seconds: values will be deleted from the selected User memory.

You can delete data from both User Memories by removing the batteries from the monitor.

## Setting the date and time

The date and time can be set using the two buttons on the top of the monitor.

1. press the  button : a flashing digit corresponding to the month will appear on the display.
2. use the  button to set the month; the month will change each time you press the button. After selecting the month, press the  button to confirm. A flashing digit corresponding to the day will now appear on the display.
3. set the day, hour and minutes, following the procedure in step 2. Use the  button to edit your settings and the  button to confirm them.
4. after setting the date and time, "0" will appear on the display. The monitor is ready for use.

## Troubleshooting

If any faults should occur during use, please check the following:

Symbol on the display	Likely cause	Remedy
The display is blank when you press the ON/OFF button	Low batteries	Replace with 2 new alkaline batteries.
	Batteries have not been inserted the right way round	Put the batteries back in, the right way round.
The EE symbol is displayed or the pressure value seems to be too low/high	Cuff fitted and/or positioned incorrectly	Put the cuff back on properly
	The user talked or moved while taking the measurement	Measure blood pressure again, without moving.
	The user moved his/her wrist during the measurement.	
The symbol E1 is displayed	Problem with the air circuit	Repeat the measurement
The symbol E2 is displayed	Pressure above 300 mmHg	Turn off the monitor, turn it back on and repeat the measurement
The symbol E3 is displayed	Data error	Remove the batteries, put them back in and repeat the measurement

## Warning:

1. the unit contains high-precision components. Therefore avoid extreme temperatures, humidity and direct sunlight. Avoid dropping or vigorously shaking the main unit and protect it from dust.
2. clean the monitor body and cuff carefully with a soft, slightly damp cloth. Do not press. Do not fold the pre-shaped cuff and do not wash with chemical detergents or use thinners, medical alcohol or petrol as cleaning fluids.
3. remove the batteries if the unit is not used for long periods.
4. do not press the ON/OFF button unless the monitor is on the wrist.
5. do not let children use the monitor.
6. do not disassemble the unit or cuff.
7. if the unit is stored at temperatures near freezing point, allow it to acclimatize at room temperature before use.
8. do not under any circumstances attempt to repair the device by yourself if it does not work. Do not try to open or disassemble the unit, or adjust anything inside it. Contact a qualified technician.
9. measurements on people with common arrhythmia (atrial or ventricular premature beats or atrial fibril-

lation), circulatory disorders, kidney problems or who have suffered a stroke or who are unconscious may not be accurate, because of their condition. In these cases, consult a doctor.

10. stored data will be deleted if the batteries are removed, so if you are monitoring your blood pressure, please note down values before removing the batteries.

### Specifications\*:

Measurement method	oscillometric method
Measuring range:	blood pressure 40 ~ 250mmHg pulse rate 40 ~ 199 beats/minute
Pressure sensor	semiconductor
Accuracy	blood pressure $\pm 3$ mmHg pulse rate: $\pm 5\%$ of reading
Inflation	pump driven
Deflation	automatic pressure release valve
Memory capacity	60 measurements from 2 users (total 120 measurements)
Auto shut-off	60 seconds after the last button is pressed
Operating environment	temperature: 10°C ~ 40°C (50°F ~ 104°F) humidity: 40 ~ 85% RH max.
Storage	temperature: -10°C ~ 60°C (14°F ~ 140°F) humidity: 10 ~ 90% RH max.
Power supply	2 AAA LR03 batteries, 1.5V DC 3V 1A DC
Size in mm	69 (L) x 69 (W) x 79 (H)
Weight	123 g (without the batteries)
Wrist circumference	13.5 ~ 22 cm (5.3" ~ 8.7")

BF-type applied part

\* These specifications are subject to change without notice.

Note: Batteries are hazardous waste — dispose of them in recycling bins. Do not dispose of them with household waste.

**This product complies with EU Directive 2002/96/EC.**

The crossed bin symbol on the appliance indicates that the product, at the end of its life, must be disposed of separately from domestic waste, either by taking it to a separate waste disposal site for electric and electronic appliances or by returning it to your dealer when you buy another similar appliance. The user is responsible for taking the appliance to a special waste disposal site at the end of its life. If the disused appliance is collected correctly as separate waste, it can be recycled, treated and disposed of ecologically; this avoids a negative impact on both the environment and health, and contributes towards the recycling of the product's materials. For further information regarding the waste disposal services available, contact your local waste disposal agency or the shop where you bought the appliance.

## Notes

1. The device is manufactured according to applicable standards and complies with all requirements on protection relating to electromagnetic compatibility; special precautions should however be taken for EMC and the device should be installed and operated in compliance with the information on EMC in this leaflet.
2. Portable and mobile radio communication equipment may affect operation of the electromedical device. When electromagnetic disturbance occurs, repeat arterial blood pressure readings if values are uncertain.

**TABLE 201 – electromagnetic emissions**

<b>Guidance and Manufacturer's Declaration – Electromagnetic Emissions</b>		
The CS 530 is intended for use in the electromagnetic environment specified below. The customer or user of the CS 530 should ensure the device is used in this environment.		
Emissions test	Complies	Electromagnetic environment - guidance
RF emissions CISPR 11	Group 1	The CS 530 monitor only uses RF energy for its internal function, therefore its RF emissions are very low and are not likely to cause any interference in nearby electronic equipment.
RF emissions CISPR 11	Class B	The CS 530 monitor is suitable for use in all environments including domestic ones and those directly connected to the public power supply network that supplies buildings used for domestic purposes.

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