CS 510

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AUTOMATIC WRIST BLOOD PRESSURE MONITOR

Features:

- 90 automatic memories
- One Touch easy measurement
- Easy-to-read display
- Auto shut-off

Introduction

Congratulations on your purchase of the new CS 510. This fully automatic wrist blood pressure monitor is easy to use and ideal for daily measurements. The large display clearly shows systolic and diastolic pressure and pulse rate after each reading.

In addition, the built-in memory of the CS 510 enables it to store up to 90 measured values, an ideal feature for users who wish to monitor and track their blood pressure readings on a regular basis. The CS 510 is compact and portable, making it ideal for home use and travel. Blood pressure measurements taken with the CS 510 are equivalent to those obtained by a trained person using the armband/stethoscope auscultation method, within the limits prescribed by the American National Standard for Electronic or Automated Sphygmomanometers. The CS 510 is designed for home use by adults.

Attention:


Please read this manual carefully before use. For specific information on your own blood pressure, please contact your doctor. Please be sure to keep this manual.

Quick Reference

A. Install batteries

Type of batteries
2 alkaline, size
AAA

D. Press the START button

Remain still during measurement. The slightest movement could result in an inaccurate measurement.

B. Apply armband

Apply armband to left wrist with the palm facing up.


E. Automatic memory

Your new CS 510 can store up to 90 measured values.

C. Correct position

Rest your wrist on a surface, level with your heart.

F. Recall measurements stored in the memory

Press the Memory button to display previous measurement readings.

Preliminary remarks

This blood pressure monitor complies with the European regulations and bears the EC mark "0068". The quality of the device has been verified and conforms to the provisions of European Council directive 93/42/EEC dated 14th June 1993 regarding medical devices, as well as EMC directive 89/336/EEC:

EN 1060-1 (Nov. 1995) Non-invasive blood pressure measuring equipment. General requirements

EN 1060-3 (Feb. 1997) Non-invasive blood pressure measuring equipment. Supplementary requirements for electromechanical blood pressure measuring systems

To ensure accurate measurements, the calibration of the device should be checked every two years.

How the CS 510 works

The CS 510 uses the oscillometric method to measure arterial blood pressure. Before the armband is inflated, the device will establish a baseline armband pressure equivalent to the air pressure. At the press of a button, the armband is automatically inflated. As it inflates, the CS 510 will measure your arterial blood pressure, analysing its oscillations. Once the amplitude and gradient of the pressure oscillations have been determined, the device will display the systolic and diastolic pressure and pulse rate.

Arterial blood pressure

What is arterial blood pressure?

Arterial blood pressure is the pressure exerted on the artery as blood flows through it. The pressure measured when the heart contracts and pumps out blood is the systolic pressure (highest). The pressure measured when the heart dilates and blood flows back is called diastolic (lowest) blood pressure.


Why Measure Your Blood Pressure?

Amongst the various health problems afflicting modern man, problems linked with high blood pressure are by far the most common. The dangerously close relationship between high blood pressure and heart disease, characterised by high morbidity, has made it necessary to monitor blood pressure in order to identify those at risk.

Blood pressure standard

The World Health Organization (WHO) has established a blood pressure standard, identifying risk areas for hypotension (low blood pressure) and hypertension (high blood pressure). This standard, however, is a general guideline as individual blood pressure measurements vary according to population, age group, etc.


It is important to consult your doctor regularly. Your doctor will be able to tell you your normal blood pressure range as well as the point at which you should consider yourself at risk.

Blood Pressure Fluctuation


Blood pressure fluctuates all the time!

Do not worry if two or three measurements are high. Blood pressure fluctuates during the month and even throughout the day. It is also influenced by season and temperature


Carry Case:


Please store the blood pressure monitor in the carry case as shown below.

Name / function of the various parts:

Key to display Display:

Symbols

mem.

Showing when the measurement has been memorised or has been recalled from the memory.

Showing the pulse rate per minute

Showing when the batteries need replacing

Showing if a mistake is made during measurement

Power supply

Inserting batteries

1. Remove the battery cover
2. Insert the first battery in the lower position on top of the ribbon.
3. Insert the other battery and fold the ribbon over it.
4. Replace the cover by placing the end tabs of the battery cover in the slots on the top of the battery compartment and pushing down the other end to secure the battery cover.

Replace the batteries when:

1. The flat battery symbol appears on the display
2. Nothing appears on the display when the ON/OFF/START button is pressed

Note: Batteries are hazardous waste. Do not dispose of them together with domestic refuse.

Applying the armband

1. Remove all watches, jewellery, etc. before applying the wrist monitor. Roll up your sleeve, as the armband must be applied to bare skin for correct measurement.


2. Apply armband to left wrist with the palm facing up, as shown in Fig. A.

3. Make sure that the edge of the armband is positioned about 1 cm from your palm, as shown in Fig. B.

4. In order to ensure accurate measurement, fasten the Velcro strap securely around your wrist so that there is no extra space between the armband and your wrist, as shown in Fig. C.

If the armband is not tight enough, the measurement values may be distorted.

5. If your doctor has diagnosed poor circulation in your left arm, carefully place the armband around your right wrist, as shown in Fig. D.

Correct Measuring Position

1. Place your elbow on a table so that the armband is at the same level as your heart, as shown in Fig. A.


Note: Your heart is located slightly below your armpit. Relax your entire body, especially the area between your elbow and fingers.

2. If the armband is not at the same level as your heart or if you are unable to keep your arm completely still whilst the measurement is being taken, use a soft object such as a folded towel to support your arm, as shown in Fig. B. Do not allow hard objects to come into contact with the wrist armband.

3. Turn your palm upwards

4. Sit upright and take 5-6 deep breaths.

Avoid leaning back whilst the measurement is being taken, as shown in Fig. C.



Helpful Tips

Here are a few helpful tips to help you obtain more accurate readings:

- Do not measure your blood pressure immediately after eating a large meal. To obtain more accurate readings, wait at least one hour before measuring your blood pressure.
- Do not smoke or consume alcohol before measuring your blood pressure.
- Do not measure your blood pressure when you are physically tired.
- It is important to relax during measurement. Try to rest for 15 minutes before measuring your blood pressure.
- Do not measure your blood pressure if you are tense or under stress.
- Measure your blood pressure at normal body temperature. If you are feeling hot or cold, wait a while before taking a measurement.
- If the monitor is stored at very low temperature (near freezing), put it in a warm place for at least an hour before use.
- Leave at least 5 minutes between each measurement.

Taking a measurement

1. Place the armband on the wrist. Press the ON/START button. All the digits will light up to check the display functions. The checking procedure will take about 2 seconds.
2. After having displayed all the symbols, the display will show a flashing 0. The monitor is now ready to measure.
3. Once the monitor has inflated the armband to a pressure of approximately 180 mmHg, measurement will begin.

4. When the measurement is complete, the systolic and diastolic pressure and pulse rate values will be shown simultaneously on the display and automatically saved in the system's memory. The memory can store up to 90 measurements.

5. Press the ON/OFF button to turn the monitor off. If no button is pressed, the monitor will automatically switch off after 60 seconds.

This monitor will automatically inflate the armband to approximately 220 mmHg if the system detects that your body needs greater pressure in order to measure your blood pressure.

Note 1: This monitor automatically switches off 60 seconds after the last button pressed.

Note 2: To interrupt the measurement, simply press any button. The armband will deflate immediately after the button is pressed.



Note 3: Do not talk or move your arm or hand muscles during measurement.

Recalling values from the memory

1. To read the values stored in the memory, press the Memory button. The memory values will be shown on the LCD display. The last measurement will be shown first. Each time the Memory button is pressed, the previous memory value is displayed.

Note 1: The data in the memory will not be lost even when the Blood Pressure Monitor is switched off.

Note 2: The values stored in the memory can be deleted in 2 ways: o by removing the batteries; o by pressing the Memory button for longer than 5 seconds.

Troubleshooting

If any faults should occur during use, please check the following points:

Symptom	Check point	Correction
Nothing appears on the display when the ON/OFF/START button is pressed	Have the batteries run out? Are the batteries inserted correctly?	Replace them with 4 new alkaline batteries. Reinsert the batteries in the correct position.
The EE symbol appears on the display or the blood pressure value is excessively high or low	Is the armband correctly applied in the proper position? Did you talk or move during measurement? Did you move your wrist whilst wearing the armband?	Reapply the armband in the correct position Take the measurement again without moving.

Note: If the unit still does not work, take it to an authorised Artsana service centre. Under no circumstances should you disassemble and attempt to repair the unit yourself.

Cautionary notes:

1. The unit contains high-precision components. Therefore, avoid extreme temperatures, humidity and direct sunlight. Avoid dropping or vigorously shaking the main unit and protect it from dust.
2. Clean the blood pressure monitor body and armband carefully with a slightly damp, soft cloth. Do not press. Do not wash the armband or use chemical cleaners on it. Never use thinner, alcohol or petrol as cleaning fluids.
3. Leaky batteries can damage the unit. Remove the batteries if the unit is not used for long periods.
4. Do not press the ON/START button when the armband is not placed around the wrist.
5. Do not let children use the device.
6. Do not disassemble the main unit or armband.
7. If the unit is stored at temperatures near freezing point, allow it to acclimatize at room temperature before use.
8. The CS 510 is not field serviceable. Do not attempt to open the device using tools or try to adjust anything inside it. If you have any problems, please contact an authorised Artana service centre.
9. The device is not suitable for use by patients with common arrhythmia (atrial or ventricular premature beats or atrial fibrillation), diabetes, circulatory disorders, kidney problems, those who have suffered a stroke or are in a state of unconsciousness.

Specifications*:

Measurement method:	oscillometric
Measurement range:	pressure: 20-280 mmHg; pulse 40-180 beats per minute
Pressure sensor:	Semiconductor
Accuracy:	pressure: ± 3 mmHg; pulse: ± 5 % of reading
Inflation:	pump driven
Deflation:	automatic pressure release valve
Memory capacity:	90
Auto shut-off:	1 minute after last button was pressed
Operating environment:	temperature: 5°C - 45 °C (41°F - 113 °F); relative humidity: 40 - 85 %
Storage environment:	temperature: -10 °C – 60 °C (14°F - 140°F); relative humidity: 10 – 95%
DC power source:	2 x LR6 (AAA) batteries 1.5V DC 6V 1A DC
Dimensions:	79(L) x 72(W) x 71(H) mm
Weight:	125 g (gross weight) (without batteries)
Wrist circumference:	13.5 - 22 cm (5.3" – 8.7")
Max. usage time:	max. 2 min.
	BF type classification

* Specifications are subject to change without notice.

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