IT

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EN

INSTRUCTION FOR USE

Disposable sterile pen needle for subcutaneous drug injections. The equipment is meant for use by trained users or health care personnel.

WARNING: People with reduced mental/physical capacities and children should be assisted while using the device. It is advisable to always consult your doctor about the length of the needle and the most suitable injection technique to be used (angle type, with or without folding). Here are the steps to follow to ensure accurate use of the pen needle:

- Check the integrity of the seal and of the pen needle. Do not use if damaged. Check the expiry date of the device. Do not use if expired.
- Check the compatibility between the pen needle and the pen injector device to be used for administering the medication, considering the information in this manual and on the packaging.
- Remove the seal from the pen needle. (See Fig. 1)
- Connect the pen injector device through pushing the pen needle straight on the pen injector device and twist until it is tight (see Fig. 2a).
- WARNING:** Pay attention during the connection between the pen needle and the pen injector device. If the needle and the pen injector device are not properly aligned, the pre-filled drug cartridge side needle tip can bend and the medication or dose may not be injected, which may result in serious injury or death.
- Remove the plastic detachable outer cover of the needle. (See Figures A and 2b). Do not discard, you will need it later to remove the pen needle from the pen injector device.
- Remove the inner protective cap by pulling it straight out, in order to prevent bending of the needle and/or damaging of the needle tip. (See Figures B and 3).
- Always keep a spare pen needle as a replacement.

WARNING: Remove both the outer cover and the inner needle cover before an injection (See Fig. A and Fig. B, respectively). If both the outer cover and the inner needle cover are not removed before use, the medication or dose may not be injected, which may result in serious injury or death (See Fig. 4).

- Prime your pen injector device as per your pen's instructions for use, so as to prevent incorrect dosage caused either by air in the cartridge or a clogged needle. If no drug flow is seen after several attempts, attach a new pen needle. Always keep a spare pen needle as a replacement.
- Set your drug dose.
- Make the injection slowly according to the technique and instructions suggested by your healthcare professional (see Fig. 10) performing the injection in a different location each time (see Fig. 9), in order to avoid lipodystrophy. Do not perform the injection in site affected by lipodystrophy. Refer to the instruction on the drugs fU (e.g. drugs temperature). Do not withdraw the needle immediately after the injection, keep it at least for 10 seconds. (See Fig. 6) If after the needle is removed from the skin, you notice medication dripping, keep the needle inserted longer the next time, at the same dosage. Note: do not massage the injection site after the dose delivery.

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Fig. 10

WARNING: Do not change the direction of the pen needle while it remains in the body, as this can result in bending or breaking of the needle.

• At the end of the injection, the outer cover may be used to cover the needle to avoid accidental punctures during the needle unscrewing and its disposal. (see Fig. 7)

WARNING: attention should be taken during covering the needles. For needles longer than 6 mm, like G29x12 mm, do not re-cap, in order to avoid that cannula gets stuck in the outer cover.

• Unscrew the pen needle from the pen injector device after each injection, in order to prevent air from entering the cartridge and the drug from leaking out (see Fig. 8). Finally dispose it in a sharp container respecting applicable regulations and laws.

• This device is single use, single patient only.

WARNING: in case of improper uses:

- If the needle is left or if the needle has been reused on the pen between injections, it may cause:
 - Breakage of the needle.
 - Inaccurate dosage since the composition of the medication may be altered.
 - Air getting into the cartridge through the needle can reduce the normal flow of medication, causing a reduced drug's delivery.
- The drug's crystallization/chemical degradation, with a partial or complete needle's obstruction, causing difficulty in or inability to carry out the injection.
- Contracting infections due to loss of sterility.
- Damage of the needle's tip ("hook" shaped) and loss of needle lubrication, causing an increased pain and possible tearing of the skin.

WARNING: Report any serious accident occurring in relation to the device to your competent authority and to the manufacturer.

Store in a fresh and dry place, far from chemical substances.

Manufacturing date: first 4 digits of LOT (e.g. 23039999; 23=2023; 03=March)

LOT Batch code

Single use. Do not re-use

REF Catalog number

Caution

Use-by date

Do not use if package is damaged and consult instructions for use

MD Medical device

Sterilized using ethylene oxide
Single sterile barrier system

Cannot be re-sterilized

Packaging material identification for recycling

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| | |
|-----------------|--|
| AstraZeneca | Byetta Pen 5 mcg, Byetta Pen 10 mcg, Symlin Pen 60, Symlin Pen 120 |
| Berlin Chemie | BerlinPen Areo 2, BerlinPen Areo 3 |
| BD | Wangbangpen |
| Bioton | GensuPen, GensuPen Improve, Genupen 2 |
| Brighter | Actiste |
| Eli Lilly | HumaPen Luxura, HumaPen Luxura HD, HumaPen Memoir, Humalog KwikPen, Humalog Mix 50/50 KwikPen, Humalog Mix 75/25 KwikPen, Humalog U200 KwikPen, Basaglar KwikPen, HumaPen Savvio, Forsteo, Humulin N KwikPen, Humulin R KwikPen, Humulin 70/30 KwikPen |
| Owen Mumford | AutoPen Classic 1-21, AutoPen Classic 2-42, AutoPen 24 1-21, AutoPen 24 2-42 |
| Ypsomed | Ypsopen |
| Mylan | Semglee |
| Novo Nordisk | NovoPen 3, NovoPen 4, NovoPen 5, NovoPen Echo, NovoPen Junior, Novo Rapid Flexpen, Levemir Flexpen, Novolog Flexpen, Novolog Mix 70/30 Flexpen, Innolet, Victoza, Flex Touch, FlexTouch Tresiba Pen U100, Flex Touch Tresiba Pen U200 |
| Sanofi | AllStar, AllStar Pro, ClickStar, JuniorStar, Lyxumia Adlyxin Pen, Suliqva 10-40 SoloStar, Suliqva 30-60 SoloStar, Sotliqua 100/33, Lantus SoloStar, Apidra SoloStar, Insuman Basal SoloStar, Insuman Rapid SoloStar, Insuman CombZ5 SoloStar, Toujeo SoloStar, Toujeo Max SoloStar, Toujeo DoubleStar, Insulin lispro SoloStar, Admelog SoloStar, Insulin Aspart, TactiPen |
| Tonghua Dongbao | Gansulin Pen |

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PICARE S.p.A.